# The Ahmed<sup>®</sup> Glaucoma Valve Tube Extender

Manufactured by

NEW WORLD  
MEDICAL

USA

CE  
0459

**Caution:** Federal (USA) law restricts this device to sale or on the order of a physician.

# Index

English (en).....	4
Français (fr).....	8
Deutsch (de).....	12
Italiano (it).....	16
Español (es).....	20
Português (pt).....	24

Uncontrolled Copy if Printed

Uncontrolled Copy if Printed

English (en)

# The **Ahmed® Glaucoma Valve** **Tube Extender**

## Product Information

### DESCRIPTION

The Ahmed™ Glaucoma Valve Tube Extender (AGV-TE) is an accessory designed to extend the tube of glaucoma drainage devices (GDDs). The AGV-TE allows for tube extension when the GDD tube is cut too short at the time of surgery or when postoperative complications require tube lengthening. In many situations, tube extension with the AGV-TE provides an alternative to replacing the existing GDD or placing a new GDD in a second location. The AGV-TE is designed to fit the existing tubes of glaucoma drainage devices including the AGV-FP7 and AGV-FP8 implants. The AGV-TE consists of a length of silicone tube attached to a silicone plate. The attached silicone tube has the same inner and outer diameter as the AGV implants listed above. On the opposite end of the silicone plate there is a protruding leak resistant junction to connect to the shortened tube from the GDD. The silicone plate also has two scleral fixation holes for easy suturing. The device is provided sterile and is for single use only.

### INTENDED USE

The AGV-TE is intended for glaucoma drainage tube lengthening. Tube lengthening may be required when a GDD tube is cut too short at the time of surgery or postoperatively in cases of tube retraction, inadvertent tube cut, kinked tube, tube blockage, tube dislocation, tube malpositioning, tube revision, and tube exposure which require tube extension.

### CONTRAINDICATIONS

- a. Bacterial conjunctivitis
- b. Bacterial corneal ulcers
- c. Endophthalmitis
- d. Orbital Cellulitis
- e. Bacteremia or septicemia
- f. Active scleritis
- g. No light perception

### COMPLICATIONS AND ADVERSE REACTIONS

Complications and adverse reactions during or following surgery may include:

Corneal Edema

Corneal Touch

Iris/Tube Touch

Synechiae

Exposed Scleral Graft

Choroidal Detachment

as well as known complications of aqueous shunts and general intraocular surgery including hypotony, shallow chamber, choroidal effusion, suprachoroidal hemorrhage, retinal detachment, cataract, conjunctival buttonhole, phthisis bulbi, bullous keratopathy, uveitis and endophthalmitis.

## SURGICAL PROCEDURE

Refer to the illustration of each step. The steps described below are intended as a guideline only and do not represent recommended treatment for any particular patient. The use of any specific surgical technique or maneuver is at the sole discretion of the surgeon. Surgeons should be familiar with the use of GDDs and post-operative care considerations before implanting any drainage device or accessory. Reference papers and surgical videotapes are available upon request.

1. If needed, expose the existing GDD tube from the overlying tissue and remove it from the anterior chamber (AC)
2. Examine the AGV-TE and flush it with 1 mL balanced salt solution (BSS) or sterile water using a blunt 26-gauge cannula (**figure 1**). The AGV-TE can be grasped with non-toothed conjunctival forceps.
3. Cut the existing tube to the appropriate length such that when the cut tube is inserted into the proximal end of the AGV-TE plate, the AGV-TE plate is at least 4-5 mm away from the limbus. Make sure to leave at least a 2 mm extension of the existing GDD tubing from the capsule or end of the device endplate to slide into the proximal end of the AGV-TE plate. (**figure 2**)
4. Insert the GDD tube into the proximal end of the AGV-TE plate with forceps until the inserted tube can go no farther. There is no need to suture the existing tube to the junction. (**figure 3**)

**NOTE:** An alternative technique is to slide a long 30 gauge cannula through the distal end of the AGV-TE tube, insert it into the cut tube of the GDD and then slide the AGV-TE plate down the cannula over the cut GDD tube.

**CAUTION:** Take care not to pull or dislodge the existing tube out of the proximal end of the AGV-TE plate during subsequent steps of the surgical procedure.

5. Secure the AGV-TE to the sclera with nonabsorbable sutures through the AGV-TE plate fixation holes. Rotate the suture knots into the the suture hole openings.
6. Cut the tube of the AGV-TE to the appropriate length to permit a 2-3 mm insertion of the tube into the AC. The tube should be bevel cut to an anterior angle of 30° to facilitate insertion and avoid iris incarceration. (**figure 4**)
7. Place the tube of the AGV-TE into the AC. As desired, enter the AC through an existing needle track or a newly created needle track. when creating a new entry pathway into the AC, enter 1-2mm from the limbus with a sharp 23-gauge needle to create a needle track, parallel to the iris. Some surgeons prefer to enter the AC from at least 3mm away from the limbus.


**NOTE:** Using an existing track may increase the chances of early leakage around the tube. Also, watertight suturing of the opening of the existing track will help prevent unwanted leakage of aqueous humor.

**CAUTION:** Care must be taken to ensure that the tube of the AGV-TE does not contact the iris or corneal endothelium after insertion in the event a new needle track is created.


8. Cover the exposed AGV-TE plate and extension tube with a piece of preserved donor sclera, pericardium, cornea, or other suitable patch graft material that is sutured into place.
9. Close the conjunctiva

**NOTE:** As an alternative to Step 6, a 2/3 thickness limbal-based scleral flap may be made. The tube is inserted into the AC through a 23-gauge needle puncture made under the flap. The flap is sutured closed.


**Figure 1**


**Figure 2**

**Figure 3**

**Figure 4**


### SYMBOLS USED ON PACKAGING

Symbol	English	Symbol	English
	Manufacturer		Sterilized using irradiation
	Authorised representative in the European Community		Batch Code
	Do not re-use		Do not re-sterilize
	Date of Manufacture YYYY-MM		Serial Number
	Consult Instructions for Use		Do not use if package is damaged
	Use-by date		


## AGV-TE

Trade Name: Ahmed® Glaucoma Valve Tube Extender

### Top View

### Side View

**Plate:** Medical Grade Silicone

**Extension Tube:** Medical Grade Silicone

**Max Plate Thickness:** 1.1 mm

**Max Plate Length:** 4.8 mm

**Max Plate Width:** 3.0 mm

**Junction Inner Diameter:** 0.635 mm

**Extension Tube Length:** 25.4 mm

**Extension Tube Outer Diameter:** 0.635 mm

**Extension Tube Inner Diameter:** 0.305 mm

## WARNINGS, PRECAUTIONS

Do not use the device if sterile package integrity has been compromised. Do not re-sterilize the implant. Before using the AGV-TE, the implanting surgeon should be skilled in glaucoma filtering procedures and familiar with the use of drainage devices, as well as post-operative care required.

Flush the AGV-TE with fluid as described in the instructions to purge any air from the lumen of the device.

Using a needle to purge the AGV-TE of air can puncture the tube resulting in an undesirable leak or inability to prime. Only use a blunt cannula to purge the AGV-TE of air.

The Ahmed™ Glaucoma Valve Tube Extender is intended for single use, and should not be reused. If the AGV-TE is explanted or not used once removed from sterile packaging, it shall be discarded to avoid deterioration, cross contamination and/or infection.

## HOW SUPPLIED

The implant is supplied sterile in a sealed pouch. Product Information, Patient Chart Labels, Implant Notification Card, and a Patient ID Card are also enclosed with the sterile package. The implant has been terminally sterilized by Gamma Radiation. Sterility is assured provided that the peel pouch has not been opened or damaged, and the sterility expiration date has not lapsed. The peel pouch is intended to be opened using sterile technique, allowing the implant to be dropped into the sterile field.

NOTE: The manufacturer disclaims all warranties expressed or implied, including but not limited to suitability for a particular purpose.

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