Why is this information necessary?

This patient information leaflet serves to inform patients on what to expect regarding the required care, precautions, and follow-up before and after Ahmed ClearPath® surgery.

What I need to know before glaucoma treatment surgery

What is glaucoma?

Glaucoma is an eye disease that damages the optic nerve of the eye. If glaucoma is untreated, individuals are at risk of permanent vision loss or blindness. Glaucoma can happen in adults and children, and it can be present since birth or caused by different factors.

What is the Ahmed ClearPath® implant and what does it do?

Ahmed ClearPath® is a non-valved glaucoma implant designed to drain fluid buildup from high eye pressure. The tube drains fluid out of the eye, thus reducing eye pressure. The implant is placed outside the eye; however, it is covered by the skin of the eye so it cannot be seen or felt.

The implant consists of a silicone tube and plate. The plate fits to the shape of the eye and provides an area where fluid can be dispersed.

What is the Ahmed ClearPath® implant used for and when is it recommended?

Ahmed ClearPath® implant is a surgical device intended to reduce eye pressure by draining fluid through a tube.

Ahmed ClearPath® implant is recommended for people with moderate to severe glaucoma and who have not responded well to other treatments. These people have usually had glaucoma surgery before, are older, and have other health problems that may affect how well the surgery works.

What I need to know after glaucoma treatment surgery

What must I generally be aware of after surgery?

Consult your physician for instructions that should be followed starting directly after your surgery. Immediately contact your physician if you experience any allergic reactions and/or adverse effects.

Ahmed ClearPath® is comprised of materials with a long history of safe clinical use. Once implanted, it remains permanently in the patient's eye and is expected to be stable over the lifetime of the patient.

It is safe to undergo Magnetic Resonance Imaging (MRI) after treatment with the Ahmed ClearPath® implant.

What warnings and precautions should I consider?

Here are some warnings and precautions to consider for ensuring a safe and successful recovery:

- Avoid activities that could cause strain on your eyes, such as lifting heavy objects, bending, or vigorous exercise for several weeks after surgery.
- Avoid rubbing or touching your eyes, as this can damage the implant or increase the risk of infection.
- Use prescribed eye drops as directed by your surgeon to prevent infection and inflammation.
- Attend all follow-up appointments with your eye doctor to monitor your progress and make any necessary adjustments to your treatment plan.
- Be aware of potential signs of infection or complications, such as pain, redness, fever, eye swelling, or vision changes. Contact your doctor immediately if you experience any of these symptoms.

- Avoid swimming, hot tubs, or any other activity that could expose your eyes to water for several weeks after surgery to prevent infection.
- Wear protective eyewear, such as goggles or sunglasses, to shield your eyes from dust, wind, or other irritants.
- Notify your doctor if you experience any sudden changes in vision like cloudiness or severe eyepain, as these may be signs of an implant malfunction or other complications.

Remember, it is essential to follow your doctor's instructions carefully and communicate any concerns or questions you may have about your recovery.

What are the possible risks and potential post-surgery adverse reactions?

Based on the health history of the patient, and the nature of glaucoma drainage device procedures, complications and adverse reactions may occur. After the surgery, your surgeon will examine your eye and will evaluate if any of these potential events or risks are present. Glaucoma surgeons are well-trained to treat the various complications that can occur after glaucoma surgery.

While surgery is generally safe, there are some general possible complications that can occur. These include but are not limited to:

- Bleeding
- Tube obstruction
- Device migration
- Corneal damage
- Swelling of the iris

As well as known complications of other types of glaucoma surgery and general eye surgery, such as:

- Low pressure inside the eye
- Serious infection inside the eye
- Clouding of the lens inside the eye (cataract)
- Detachment of the retina
- Inflammation inside the eye

It is important to note that these complications are rare, and most people who undergo glaucoma implant surgery experience a significant improvement in their condition. However, it is important to be aware of the risks and to discuss them with your doctor.

Overall benefit-risk evaluation

There is no known cure for glaucoma, and individuals with untreated glaucoma are at risk of permanent vision loss or blindness. Lowering eye pressure is the only adjustable risk factor for minimizing glaucoma progression.

The overall goal of the Ahmed ClearPath® is the immediate reduction in eye pressure and/or reducing the need for pressure-lowering medications. This clinical benefit may reduce over time due to the natural wound healing around the implant which may restrict fluid to the surrounding tissue.

Most complications and adverse reactions after surgery are not intrinsic to the implant itself but are due to the nature of the surgery and the presence of a foreign body in the eye.

The benefit of preserving vision by using the Ahmed ClearPath® outweighs the potential complications and adverse reactions, since progression of glaucoma can lead to significant loss of vision and permanent blindness.

There is no single target eye pressure that is appropriate for every patient; target pressure needs to be estimated separately for each eye and adapted to changing patient conditions.

Talk to your doctor to establish the eye pressure target suitable for you; this target pressure may be updated at each monitoring visit based on changes in your medical condition. New World Medical has considered and mitigated as far as possible all significant risks and deemed them acceptable.

Further Information:

Contents: The implant consists of a medical-grade silicone tube and plate. The silicone plate is barium filled and contains a minimum of 8% and a maximum of 10% of barium sulfate (BaSO_4) to increase ultrasound resolution and identification with CT scans, MRIs, and plain skull films.

Model: CP250, CP350

Please refer to our Patient Implant Card for the model number of your implant.

Therapeutic Goods Administration:

Any serious incident that occurs in relation to the device should be reported directly to the local representative or the manufacturer at the email address surveillance@newworldmedical.com, as well as the Therapeutic Goods Administration at www.tga.gov.au

European Union Member States:

For customers from countries within the EU, any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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New World Medical, Inc.